

Fast, safe, non-contact measurement of human body temperature

- In Stock in the USA
- Forehead Non-Contact
- Ships the Next Day
- 1 Year Warranty, CE, RoHS

Simply point the non-contact IR (infra-red) thermometer at the subject's forehead from a distance of 2" to 6", pull the trigger and, within 1 second, obtain immediate, highly accurate readings of human body temperature. Automatically adjusts for the difference between the surface temperature of the forehead and internal body temperature of the subject, for highly accurate reading within +/- 0.3°C. Also features an audible max/min alarm, backlit display, 32 memories, and auto-power-off. Comes ready to use with one 9V battery.

Calibrated version available.



No.	Description
-----	-------------

- | | |
|---------|-----------------------------------|
| 800120 | Clinical IR Thermometer |
| 800120C | Clinical IR Thermometer Certified |

	Range	Resolution	Accuracy
Temperature °F	89 to 110.0°F	0.1°F	±0.6°F
Temperature °C	35 to 42.0°C	0.1°C	±0.3°C